# Levenhuk Rainbow 2L / D2L / 2L PLUS Microscopes


User Manual  
Návod k použití  
Bedienungsanleitung  
Guía del usuario  
Instrukcja obsługi  
Инструкция по эксплуатации

*Radost zaostřit  
Zoom ran und hab Fun!  
Amplie y disfrute  
Radość przybliżania  
Приближает с удовольствием*

**levenhuk**<sup>®</sup>  
Zoom&Joy

USA • Canada • Czech Republic • Estonia • Germany • Iceland • India • Latvia • Lithuania  
Malaysia • Netherlands • Poland • Russia • Slovakia • Spain • Turkey • Ukraine • United Kingdom

**EN**

- 1 Eyepiece
- 2 Monocular head (eyepiece tube)
- 3 Revolving nosepiece with objectives
- 4 Upper illumination
- 5 Focusing knob
- 6 A/C adapter socket
- 7 Upper illumination switch
- 8 Specimen holders
- 9 Stage
- 10 Condenser and diaphragm disk
- 11 Lower illumination
- 12 Lower illumination switch
- 13 Base

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## General use

The Levenhuk Rainbow microscope is safe for health, life and property of the consumer and the environment when properly used, and meets the requirements of international standards. The microscope is designed for observing transparent and opaque objects in the transmitted and reflected light using the bright field method, for biological use and school presentations. May be used by kids over 5 years old.

The package includes a Levenhuk K50 experiment kit that has everything you need to prepare microscope slides yourself. The detailed guide will teach you how to properly prepare specimens and conduct amazing experiments. The digital camera (for Levenhuk Rainbow D2L) was exclusively designed for use with this microscope.

### *The package includes:*

- Microscope
- Objectives: 4x, 10x and 40x (40xs for 2L PLUS model)
- Eyepiece: WF10x (WF16x for 2L PLUS model)
- Stage with specimen holders
- Diaphragm disk
- Condenser
- Built-in upper and lower LED illumination
- 3\*AA batteries
- Levenhuk K50 experiment kit
- Digital camera 0.3M (for D2L model)
- Levenhuk TouView image processing software (for D2L model)
- USB cable (for D2L model)
- User manual and warranty
- A/C Adapter 220 V, 50 Hz

### *The Levenhuk K50 experiment kit includes:*

- "Attractive Microscope. Scrutinizing the Microcosm" User Guide
- Forceps
- Hatchery for Artemia (brine shrimp)
- Microtome
- Flask with yeast
- Flask with glue
- Flask with sea salt
- Flask with brine shrimp
- Blank glass slides (5 pieces)
- 5 samples (fly's limb, onion, cotton stem, slice of a tree trunk, pine)
- Pipette
- Dust cover

## Microscope parts

**Base.** It supports the weight of the microscope and houses the illumination source, electronics and control mechanisms.

**Eyepiece tube.** Combines the eyepiece with the objectives system. Holds the eyepiece, Barlow lens (below the eyepiece) or a digital camera (in place of eyepiece).

**Eyepiece and objective.** Consist of lenses that allow magnifying the image. The total magnification is calculated by multiplying the eyepiece magnification to the objective magnification.

**Revolving nosepiece.** This triple nosepiece with 3 pre-installed objectives allows you to change objectives smoothly and easily.

**Stage.** Sturdy and reliable stage with two specimen holders can be used to move your slides while observing them. The lower illumination light passes through the opening in the middle of the stage.

**Diaphragm disk.** Is located below the stage and has apertures of various diameter to adjust the passing light rays. Rotate the disk to select the desired aperture.

**Condenser.** A system of lenses located below the stage that allows gathering and directing of the light rays.

**Focusing knob.** A coarse focusing system allows moving the stage up and down adjusting the sharpness of the specimen image.

**Upper and lower illumination.** LED illumination with adjustable brightness. The illumination system can be powered by AA batteries or by A/C. The upper illumination is used to observe opaque objects whereas the lower one allows observing transparent objects. Use the both illuminations to study semi-transparent objects.

**A/C adapter socket.** Is used to connect the microscope to A/C power.

## Using the microscope

### Getting started

- Unpack the microscope and make sure all parts are available.
- Move the stage to the lowermost position using the focusing knob.
- Make sure the batteries are correctly installed in the battery compartment; insert new batteries if required. Or, you can connect the microscope to A/C power using the A/C adapter.
- Insert the eyepiece into the eyepiece tube.
- Slowly adjust the illumination brightness, from dark to light.

### Focusing

- Place a specimen on the stage and fix it with the holders.
- Select the 4x objective rotating the revolving nosepiece.
- Move the specimen to place its thickest part exactly under the objective.
- Rotate the focusing knob to slowly raise the stage until the objective is close to the specimen; keep checking the distance between the objective and the object to avoid their contact. **CAUTION:** The objective should not touch the specimen, otherwise the objective or/and the specimen might be damaged.
- Look through the installed eyepiece and lower the stage slowly rotating the focusing knob until you see the specimen image.
- Such adjustment protects the frontal lens from contacting the object when you use objectives of other magnifications; though, slight refocusing might be required.
- If the image is too bright, rotate the diaphragm disk until the passing light ray is reduced to a comfortable brightness level. If the image is too dark, select a larger aperture to increase the light ray.

### Selecting the objective

Start your observations with the lowest magnification objective and select a specimen segment for detailed research. Then move the specimen to center the selected segment in the field of view, to make sure it keeps centered when the objective is changed to a more powerful one. Once the segment is selected, you should center its image in the microscope's field of view as precisely as possible. Otherwise, the desired segment might fail to center in the field of view of the higher power objective. Now you can switch to a more powerful objective by rotating the revolving nosepiece. Adjust the image focus if required.

## Digital camera

The digital camera is installed in the eyepiece tube in place of the eyepiece. It allows to observe specimens in fine detail and true colors on your PC monitor and save images on the hard drive. Levenhuk ToupView software allows viewing and editing the object images.

Megapixels	0.3
Video recording option	yes
Mounting location	eyepiece tube (replaces the eyepiece)
Image format	*.bmp, *.jpg, *.jpeg, *.png, *.tif, *.tiff, *.gif, *.psd, *.ico, *.emf, etc.
Exposure	ERS
White balance	auto/manual
Exposure control	auto/manual
Software	Levenhuk ToupView
Software image control	image size, brightness, exposure control
Port	USB 2.0, 480Mb/s
System requirements	Windows XP/Vista/7/8/10 (32 and 64 bit), processor up to 2.8GHz Intel Core 2 or higher, USB 2.0 port
Camera power supply	USB cable

## Specifications

	2L	2DL	2L PLUS
Type	biological	digital	biological
Magnification, x	40–400	40–400	64–640
Head	monocular, rotatable 360°, inclined 45°		
Optics material	optical glass		
Body material	plastic	plastic	metal
Eyepiece tube diameter, mm	23.2		
Eyepieces	WF10x	WF10x	WF16x
Revolving nosepiece	3 objectives		
Objectives	4x, 10x, 40x	4x, 10x, 40x	4x, 10x, 40xs (spring)
Stage, mm	90x90, with specimen holders		
Stage moving rage, mm	0–15, vertical	0–15, vertical	0–8, vertical
Focusing	coarse		
Condenser	NA 0.65, diaphragm disk (6 apertures)		
Illumination	upper and lower LED illumination with adjustable brightness		
Power source	220V 50Hz, or 3*AA batteries		
Levenhuk K50 experiment kit	yes		
Digital camera	no	yes	no

Levenhuk reserves the right to modify or discontinue any product without prior notice.

### Notes:

1. The Levenhuk K50 experiment kit is supplied without an individual plastic case.
2. Batteries might be pre-installed in the battery compartment by the manufacturer.

### Care and maintenance

- **Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS.**
- Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions.
- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center.
- Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws.
- Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.
- Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures.

- When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may heat up during use, and might become a fire hazard.
- Always unplug the microscope from a power source before opening the base or changing the illumination lamp. Regardless of the lamp type (halogen or incandescent), give it some time to cool down before trying to change it, and always change it to a lamp of the same type.
- Always use the power supply with the proper voltage, i.e. indicated in the specifications of your new microscope. Plugging the instrument into a different power outlet may damage the electric circuitry of the microscope, burn out the lamp, or even cause a short circuit.
- **Seek medical advice immediately if a small part or a battery is swallowed.**

### **Battery safety instructions**

- Always purchase the correct size and grade of battery most suitable for the intended use.
- Always replace the whole set of batteries at one time; taking care not to mix old and new ones, or batteries of different types.
- Clean the battery contacts and also those of the device prior to battery installation.
- Make sure the batteries are installed correctly with regard to polarity (+ and -).
- Remove batteries from equipment that is not to be used for an extended period of time.
- Remove used batteries promptly.
- Never attempt to recharge primary batteries as this may cause leakage, fire, or explosion.
- Never short-circuit batteries as this may lead to high temperatures, leakage, or explosion.
- Never heat batteries in order to revive them.
- Remember to switch off devices after use.
- Keep batteries out of the reach of children, to avoid risk of ingestion, suffocation, or poisoning.
- Utilize used batteries as prescribed by your country laws.

### **Levenhuk International Lifetime Warranty**

All Levenhuk telescopes, microscopes, binoculars and other optical products, except for accessories, carry a **lifetime warranty** against defects in materials and workmanship. Lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from date of retail purchase. Levenhuk will repair or replace such product or part thereof which, upon inspection by Levenhuk, is found to be defective in materials or workmanship. As a condition to the obligation of Levenhuk to repair or replace such product, the product must be returned to Levenhuk together with proof of purchase satisfactory to Levenhuk.

This warranty does not cover consumable parts, such as bulbs (electrical, LED, halogen, energy-saving and other types of lamps), batteries (rechargeable and non-rechargeable), electrical consumables etc.

For further details, please visit our web site: <https://www.levenhuk.com/warranty>. If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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## Levenhuk prepared slides sets

The Levenhuk slides sets include thoroughly prepared specimens of various biomaterials as well as blank slides and cover glasses for you to make your own specimens. Prepared slides, blank slides and cover slips can be used with any microscope model.

### Levenhuk Worldwide:

USA: [www.levenhuk.com](http://www.levenhuk.com)  
Canada: [ca.levenhuk.com](http://ca.levenhuk.com)  
Czech Rep.: [www.levenhuk.cz](http://www.levenhuk.cz)  
Estonia: [ee.levenhuk.com](http://ee.levenhuk.com)  
Germany: [de.levenhuk.com](http://de.levenhuk.com)  
Iceland: [is.levenhuk.com](http://is.levenhuk.com)

India: [in.levenhuk.com](http://in.levenhuk.com)  
Latvia: [lv.levenhuk.com](http://lv.levenhuk.com)  
Lithuania: [lt.levenhuk.com](http://lt.levenhuk.com)  
Netherlands: [nl.levenhuk.com](http://nl.levenhuk.com)  
Poland: [pl.levenhuk.com](http://pl.levenhuk.com)  
Russia: [www.levenhuk.ru](http://www.levenhuk.ru)

Slovakia: [sk.levenhuk.com](http://sk.levenhuk.com)  
Spain: [es.levenhuk.com](http://es.levenhuk.com)  
Turkey: [tr.levenhuk.com](http://tr.levenhuk.com)  
Ukraine: [www.levenhuk.ua](http://www.levenhuk.ua)  
United Kingdom: [uk.levenhuk.com](http://uk.levenhuk.com)  
EU: [eu.levenhuk.com](http://eu.levenhuk.com)

Levenhuk, Inc. 924-D East 124th Ave. Tampa, FL 33612 USA  
Levenhuk® is a registered trademark of Levenhuk, Inc.  
© 2006–2016 Levenhuk, Inc. All rights reserved.  
20151229

**levenhuk**  
Zoom&Joy